

| Product | Code | Packing |
|---|--|----------------------------------|
| HiShield™ - 5 in 500ml Bottle w/o pump (Refill pack) | C0135-1x10NO C0135-1x20NO | 1x10 no 1x20 no |

HiShield™ HC

HiShield™ HC is disinfectant antiseptic solution containing Chlorhexidine Gluconate Solution IP, Cetrimide IP is suitable for disinfection of skin, hard surface and sterilization of medical equipments. Superior microbicidal action against highly contagious microbes including antibiotic resistant bacteria. Effective against MRSA stains.

Microbicidal Properties: Bactericidal, Fungicidal, Effective against Yeast, Mold and Protozoa

| | | |
|--|--------------------------------|----------------|
| HiShield™ HC In 500ml bottle | New C0191-1x20NO | 1x20 no |
| HiShield™ HC In 1 litre can pack | New C0192-1x6NO | 1x6 no |

| Product | Code | Packing |
|-------------------------------|------|---------|
| HiShield™ PVP Solution | | |

HiShield™ PVP Solution is a topical Microbicide containing Povidone Iodine IP : 10% W/V (1% W/V available iodine). Its microbicidal activity is retained in presence of blood pus & serums.



| | | |
|--|--|----------------------------------|
| HiShield™ PVP Solution In 500ml Bottle | C0144-1x10NO C0144-1x20NO | 1x10 no 1x20 no |
| HiShield™ PVP Solution In 1 lit Can Pack | C0145-1x5NO | 1x5 no |
| HiShield™ PVP Solution In 5 lit Can Pack | C0146-1x2NO | 1x2 no |

Instruments / Equipments (Disinfection & Cleaning)

AlkaSept™

AlkaSept™ is a residue-free, exceptionally material-compatible disinfectant for chemo-thermal reprocessing of endoscopes. The product achieves a high micro-organism reduction and offers staff and patients an effective protection against infection. It contains Active ingredient as Glutaraldehyde & other ingredients: complexing agents, corrosion inhibitors and activator.



| | | |
|---------------------------------------|--------------------|---------------|
| AlkaSept™ In 1 lit Can Pack | C0152-1x5NO | 1x5 no |
| AlkaSept™ In 5 lit Can Pack | C0102-1x2NO | 1x2 no |

AlkaSept™ Plus

AlkaSept Plus™ is a formaldehyde free surface disinfectant cleaner with active ingredients Benzalkonium chloride IP, Glutaraldehyde, Didecyl-dimethylammonium chloride. It is recommended for the disinfectant cleaning of medical equipment in hospitals and residential/nursing homes.



| | | |
|--|--|----------------------------------|
| AlkaSept™ Plus In 500ml Bottle | C0103-1x10NO C0103-1x20NO | 1x10 no 1x20 no |
| AlkaSept™ Plus In 1 lit Can Pack | C0104-1x5NO | 1x5 no |
| AlkaSept™ Plus In 5 lit Can Pack | C0105-1x2NO | 1x2 no |

For your requirement of **Cell Culture, Plant Tissue Culture & Molecular Biology** Laboratory Aids products, refer respective section.

| Product | Code | Packing |
|----------------------|------|---------|
| AlkaSept™ OPA | | |

AlkaSept™ OPA is a concentrated High Level Disinfectant (HLD) solution for use in re-processing heat sensitive medical devices. Its active ingredient is Ortho-phthalaldehyde along with other inert ingredients that provide disinfection in 12 minutes at room temperature for manual processing. It is effective against Bacteria (particularly Mycobacteria), Fungi and spores of Clostridium.



| | | |
|--|--|------------------------------|
| AlkaSept™ OPA In 1 lit can pack | C0161-1NO C0161-1x3NO | 1 no 1x3 no |
| AlkaSept™ OPA In 5 lit can pack | C0162-1NO | 1 no |
| AlkaSept™ OPA Solution Test Strips For potency testing of AlkaSept™ (Code no. C0161 & C0162), 50 strips per kit. | K116-1KT | 1 kt |

AlkaSept™ Active

AlkaSept™ Active is a glutaraldehydebenzalkonium concentrate suitable for cleaning surface disinfection of washable surfaces such as medical devices, functional areas of hospitals/nursing homes, laboratories and cosmetic industries using the wet-wipe procedure. It contains Glutaraldehyde, Benzalkonium chloride, 1,6 Dihydroxy 2-5 Dioxahexane and Alkyl Urea derivative.



| | | |
|---|--------------------|---------------|
| AlkaSept™ Active In 500 ml bottle | C0163-1x5NO | 1x5 no |
|---|--------------------|---------------|

| Product | Code | Packing |
|--|------------------|-------------|
| AlkaSept™ Active In 1 lit can pack | C0164-1NO | 1 no |

HiQuat Citrus

HiQuat Citrus is Disinfectant cleaner with low use concentration containing Benzalkonium chloride IP, Didecylidimethyl ammonium chloride, Citrus oil. Microbicidal Properties: Bactericidal, Fungicidal, Yeasticidal and Virucidal. Also effective against mold and protozoa.

| | | |
|---|-------------------------|--|
| HiQuat Citrus In 500ml bottle | New C0200-1x20NO | |
| HiQuat Citrus In 1 lit bottle | New C0201-1x6NO | |
| HiQuat Citrus In 5 lit can pack | New C0202-1x2NO | |

SteRapid Disinfectant Spray* for Medical Devices

Clear transparent liquid with Benzalkonium chloride IP, Isopropyl alcohol v/v as antimicrobial agents leaving no residues after application to **medical devices, surgical equipments, instruments** and other critical parts of hygienic importance. It is recommended for aerial fogging.



| | | |
|--|---|---------------------------------|
| SteRapid Spray Bottle in 500 ml Pack | C0028-1x5NO C0028-1x20NO | 1x5 no 1x20 no |
|--|---|---------------------------------|

| Product | Code | Packing |
|-------------------|------|---------|
| InstaCull™ | | |

InstaCull™ is recommended for disinfection in various areas of hospitals like NICU purpose, Respirators and AKD Units, Dental chairs, equipments, hand pieces, Air space between patients and surgeries Hospital beds, Table tops, Bed pans etc. In Pharmaceutical, Food and Dairy Industry for Laminar Air Flow, Biosafety cabinets, Table tops, Work Stations air and Surface disinfection in critical areas. It contains Isopropyl Alcohol, Benzalkonium chloride IP, Propanol, Octyldecyldimethyl ammonium chloride.



| | | |
|--|--|----------------------------------|
| InstaCull™ In 500ml Dispenser Bottle w/ pump | C0109-1x10NO C0109-1x20NO | 1x10 no 1x20 no |
| InstaCull™ In 500ml Bottle w/o pump (Refill pack) | C0110-1x10NO C0110-1x20NO | 1x10 no 1x20 no |

PowerCull™ Active

“**PowerCull Active™** is an aldehyde free instrument disinfectant for heat sensitive and heat resistant instruments with powerful cleaning effect. It contains Active ingredients as Benzalkonium chloride IP, didecyldimethyl-ammoniumchloride, along with surfactants, solvents, corrosion inhibitors, PH value regulators, Antifoaming agents.”



| | | |
|---|--|----------------------------------|
| PowerCull™ Active In 500ml Bottle | C0122-1x10NO C0122-1x20NO | 1x10 no 1x20 no |
|---|--|----------------------------------|

| Product | Code | Packing |
|------------------------|------|---------|
| PowerCull™ Plus | | |

| | | |
|---|-------------------------|---------------|
| PowerCull™ Active In 1 lit Can Pack | C0123-1x5NO | 1x5 no |
| PowerCull™ Active in 5 lit Can Pack | C0124-1x2NO | 1x2 no |
| HiShield™ HC In 500ml bottle | New C0191-1x10NO | |
| HiShield™ HC In 1 litre can pack | New C0192-1x10NO | |

PowerCull™ Plus

PowerCull™ Plus is an aldehyde free instrument disinfectant for heat sensitive and heat resistant instruments with powerful cleaning effect. Active ingredients are: Benzalkonium chloride IP, Didecyldimethyl-ammonium chloride & Poly Hexamethylene Biguanide hydrochloride.



| | | |
|---|--|----------------------------------|
| PowerCull™ Plus In 500ml Bottle | C0125-1x10NO C0125-1x20NO | 1x10 no 1x20 no |
| PowerCull™ Plus In 1 lit Can Pack | C0126-1x5NO | 1x5 no |
| PowerCull™ Plus in 5 lit Can Pack | C0150-1x2NO | 1x2 no |

HiKleen

HiKleen is a cold sterilent with excellent antimicrobial property containing per Acetic acid, Hydrogen Peroxide & Other ingredients: Acetic acid. Microbicidal properties: Bactericidal, Fungicidal, Sporicidal, Virucidal against enveloped viruses (incl HBV, HIV, HCV) Virucidal (Parvo, Adeno, Papova, Polio, Vaccinia viruses and HAV). Rapidly kills a broad spectrum of microorganisms, coliforms, yeasts and spores even at low temperature. It removes mineral deposits and organic matters.

| | | |
|-------------------------------------|------------------------|---------------|
| HiKleen In 1 lit bottle | New C0196-1x6NO | 1x6 no |
| HiKleen In 5 lit can pack | New C0197-1x2NO | 1x2 no |

For your requirement of **Cell Culture, Plant Tissue Culture & Molecular Biology** Laboratory Aids products, refer respective section.

| Product | Code | Packing |
|--------------------------|------|---------|
| HiSpark Cleaning* | | |

| | | |
|---|---------------------------------------|----------------------------|
| HiSpark Cleaning Solution, Neutral, Biodegradable Use 2-5% concentration for dilution. Special cleaner for precision instruments of glass, quartz and sensitive metals. Suitable for ultrasonic cleaning. pH: 7.2 - 7.5 | HS002-500ML HS002-5L | 500ml 5lt |
|---|---------------------------------------|----------------------------|

#: For more details refer Chemical price list

| |
|--|
| HiShield™ - AL Spray (Disinfectant) |
|--|

Sterile 70% Isopropyl alcohol solution,

HiShield™ - AL Spray, a sterile 70% IPA spray bottle of 225ml by HiMedia Labs is the ideal disinfectant for control of pathogenic and contaminating organisms in laboratories, controlled environments and process equipments.



| | | |
|---|--|----------------------------|
| HiShield™ - AL Spray Gamma irradiated, Triple pack, 5 lit Can Pack | CO154GT-1NO | 1 no |
| HiShield™ - Aerosol AL Spray Gamma irradiated, Triple pack, 225 ml in aluminium canister pack | CO158AGT-1x2NO CO158AGT-1x5NO | 2 no 5 no |
| IPA 70% In 500ml Dispenser Bottle w/Pump) | CO181-1x10NO | 1x10 no |

| Product | Code | Packing |
|--|------|---------|
| Disinfectant Swabs, Wipes & Freshner Tissue | | |

For Institutional & Industrial use. For use in Hospitals, Clinics & Offices on washable non porous surfaces of Stainless steel, Metal, Glass & Acrylic.



| | | |
|--|--|--|
| SteriSwift™ Disinfectant Wipes Size: 7" x 7" | CO100-1x2NO CO100-1x5NO CO100-1x10NO | 1x2 no 1x5 no 1x10 no |
| SteriSwift™ Disinfectant Wipes Size: 6" x 8" | CO100S-1x2NO CO100S-1x5NO CO100S-1x10NO | 1x2 no 1x5 no 1x10 no |
| HiShield™ Sterile Presaturated Wipes Size: 8" x 10" Gamma irradiated, Pouch packing Useful for clean rooms | New CO216-1x20NO | 1x20 no |
| HiShield™ Sterile Presaturated Wipes Size: 8" x 10" Gamma irradiated, Pouch packing Useful for clean rooms | New CO217-1x40NO | 1x40 no |

For your requirement of **Cell Culture, Plant Tissue Culture & Molecular Biology** Laboratory Aids products, refer respective section.

| Product | Code | Packing |
|---|-------------------------|---------|
| HiShield™ Sterile Dry Wipes Size: 8" x 10" Gamma irradiated, Pouch packing Useful for clean rooms | New C0218-1x20NO | 1x20 no |
| HiShield™ Sterile Dry Wipes Size: 8" x 10" Gamma irradiated, Pouch packing Useful for clean rooms | New C0219-1x40NO | 1x40 no |

Disinfectant Wipes in pouch packing



| | | |
|---|------------|-------|
| SteriSwift™ Disinfectant Wipes, In Pouch Packing Size: 8" x 10" | C0173-20NO | 20 no |
| SteriSwift™ Disinfectant Wipes, In Pouch Packing Size: 5" x 7" | C0174-10NO | 10 no |
| Hi-Quat Disinfectant Wipes, In Pouch Packing Size: 8" x 10" | C0175-20NO | 20 no |
| Hi-Quat Disinfectant Wipes, In Pouch Packing Size: 5" x 7" | C0182-10NO | 10 no |
| HiShield™ Disinfectant Wipes, In Pouch Packing Size: 8" x 10" | C0176-20NO | 20 no |
| HiShield™ Disinfectant Wipes, In Pouch Packing Size: 5" x 7" | C0183-10NO | 10 no |
| SteriFast™ Disinfectant Wipes, In Pouch Packing Size: 8" x 10" | C0184-20NO | 20 no |
| SteriFast™ Disinfectant Wipes, In Pouch Packing Size: 5" x 7" | C0185-10NO | 10 no |

| Product | Code | Packing |
|--------------------------|------|---------|
| Disinfectant Kits | | |



| | | |
|---|-----------|------|
| HiCare Hand & Skin Care Disinfectant Kit Kit contains: 1. HiShield™ Hand Rub - 100ml, 2. Hand Guard - 100ml, 3. HiShield™ Hand Wash - 100ml, 4. HiShield™ PVP Spray - 100ml, 5. Sterifast™ wipes (10 wipes of 5" x 7") | C0179-1NO | 1 no |
| HiCare Personal Care Disinfectant Kit Kit contains: 1. HiShield™ Hand Rub - 100ml, 2. Hand Guard - 100ml, 3. InstaCull™ Spray - 100ml, 4. HiShield™ PVP Spray - 100ml, 5. Sterifast™ Wipes (10 wipes of 5" x 7") | C0178-1NO | 1 no |
| HiCare Hand & Skin Care Disinfectant Kit NAH & SC (Non Icoholic) 1. HiShield™ NAHS - 100ml, 2. Hand Guard - 100ml, 3. HiShield™ Hand Wash - 100ml, 4. HiShield™ PVP Spray - 100ml, 5. HiShield™ wipes (10 wipes of 5" x 7") | C0180-1NO | 1 no |

Product

Code

Packing

HiPolyzyme™ Detergent Cleaner

HiPolyzyme™ Detergent Cleaner

is an enzymatic detergent cleaner for medical devices such as endoscope that is soiled with proteins, lipids, carbohydrates and starch. It contains rust inhibitor to help protect the valuable medical devices and surgical instruments. This solution is recommended to clean instruments before disinfecting it with High-level disinfectant solutions. It contains mixture of enzymes such as Protease, Lipase, Cellulase, Amylase with surfactants, Corrosion inhibitors and preservatives.



| | | |
|---|--------------------|---------------|
| HiPolyzyme™ Detergent Cleaner In 500 ml bottle | C0159-1x5N0 | 1x5 no |
| HiPolyzyme™ Detergent Cleaner In 1 lit can pack | C0160-1x2N0 | 1x2 no |

Product

Code

Packing

HiClean and Cure

HiClean and Cure is a cleaner cum stain and rust remover, containing Organic & inorganic acids, Nonionic surfactants, Corrosion inhibitors. It has no Antimicrobial property. As cleaner for entire sanitary area for rapid and reliable removal of primary minerals and organic soils. Suitable for all acid resistant material such as porcelain, Stainless Steel, chromium, Aluminum, colorfast synthetics, rubber, tiles, glass, and enamels. As anti calcification agent in the kitchen.

| | | |
|--|-------------------------|----------------|
| HiClean and Cure In 500ml bottle | New C0193-1x20N0 | 1x20 no |
| HiClean and Cure In 1 lit bottle | New C0194-1x6N0 | 1x6 no |
| HiClean and Cure In 5 lit can pack | New C0195-1x2N0 | 1x2 no |

Surface Disinfectants

HiSpark Cleaning*

Use 2-5% concentration for dilution. Universal cleaning agent for removing heavy contamination. For hard water even up to 40 d. For cleaning tables, tiles and floors in the laboratory. Tested for radioactive decontamination. Suitable for ultrasonic cleaning. pH : 11.6 - 12.0



| | | |
|--|---------------------------------------|----------------------------|
| HiSpark Alkaline Cleaning Solution, Biodegradable | HS001-500ML HS001-5L | 500ml 5lt |
|--|---------------------------------------|----------------------------|

Liquid Soap



Hiclean (liquid soap)

| | | |
|------------------------------|------------------------------------|--------------------------|
| Hiclean (liquid soap) | AS023-1L AS023-5L | 1lt 5lt |
|------------------------------|------------------------------------|--------------------------|

#: For more details refer Chemical price list